



XP-300™

Automated Hematology Analyzer

Robust and Reliable!



XP-300 – Robust and Reliable

Accurate and reliable

- Sysmex robustness for the best possible up-time
- Proven technology for accuracy of results
- Linearity to support clinical diagnostics and monitoring
- Sensitive flagging to support diagnosis by the physician

Innovative technology

- Performs rapid and accurate analysis of a 17-parameter CBC
 - Including a 3-part WBC differential with an Absolute Neutrophil Count (ANC)
 - Histograms for RBC, WBC and PLT
 - Automatic floating discriminators
- Produces accurate results comparable to other Sysmex hematology analyzers

Compact and fully integrated

- Small footprint
- Fits easily on a laboratory bench
- Ideal as a back-up for Sysmex X-Series™ and XN-Series™ analyzers

Easy operation and maintenance

- Minimal training required
- Barcoded QC and reagents
- Auto start-up and shut-down

Safe and secure

- Non-toxic, biodegradable reagent system
- Two reagents (diluent and lyse) for complete results

Network capability

- Data transfer to LIS via LAN (TCP/IP) or Serial Port (RS-232C)
- Bi-directional interface



LAB Market

لاب ماركت

- مستلزمات معملية
- كيمائيات
- أجهزة
- صيانة

01210312346
01016163354

بنها - عمارة ٢ غرب الاستناد - بجوار كنتاكي - خلف نقابة الاطباء ومطعم الرشيدى