

# DIRUI BCC-3000B Hematology Analyzer

## BCC -3000B

### Features

- 3 part differential hematology analyzer
- 18 report parameters and 3 histograms.
- Throughput: 60 samples/hour. Whole blood, pre-diluted blood.
- 1000 sample results, 570 QC results storage capacity
- Up to 14 alarm information of para-morphia cell
- Large color LED display with high resolution
- Environmental friendly cyanide-free reagents
- Weight: 42 lbs / 19 kilos
- W x D x H: 18X15X18 inches / 45X39X45 cm
- Open system: Open.



<b>Parametry mierzzone:</b>	WBC, RBC, HGB, PLT
<b>Parameters:</b>	Lym%, Mon%, Gra%, MCV, RDW-SD, RDW-CV, MPV, PDW, P-LCR, Lym#, Mon# Gra#, HCT, MCH, MCHC, PCT
<b>Histogramy :</b>	WBC, RBC, PLT
<b>Sample volume:</b>	13 uL - whole blood ; 20 uL - blood pre- diluted
<b>Efficiency:</b>	60 / h
<b>Display:</b>	5 large display , backlit LCD ( 320 x 240 )
<b>Memory:</b>	1000 results
<b>Display:</b>	5 large display , backlit LCD ( 320 x 240 )
<b>Outputs:</b>	RS-232
<b>Voltage:</b>	220 V / 110 V $\pm$ 10 % , 50 / 60Hz
<b>Dimension:</b>	410 mm x 390 mm x 460mm (HxWxD)
<b>Weight:</b>	25 kg



# LAB Market

## لاب ماركت

- مستلزمات معملية
- كيمائيات
- أجهزة
- صيانة

01210312346  
01016163354

بنها - عمارة ٢ غرب الاستاد - بجوار كنتاكي - خلف نقابة الاطباء ومطعم الرشيدى