

URIT-2900

Hematology Analyzer



| Specifications |

Parameters: WBC, LY#, MO#, GR#, LY%, MO%, GR%, RBC, HGB, HCT, MCV, MCH, MCHC, RDW-CV, RDW-SD, PLT, MPV, PDW, PCT

3-part differentiation of WBC, 19 parameters and 3 histograms

Aspiration Volume: Whole blood 18 μ L , pre-diluent 20 μ L

Throughput: Up to 30 samples per hour

Data storage: 2000 sample results with histograms

Display: 5.5 inches LCD, resolution 320X240

Dimension: 430mmX260mmX420mm **N.W.:** 14.5kg

Power requirement: AC 120 or 220V, 50/60Hz

Input/output: System keyboard, built-in thermal printer, port for an external matrix printer, RS232 to computer

Carryover: WBC, RBC, HGB, HCT, PLT \leq 2.0%

| Characteristics |

- ✔ Core parts with high quality, compact design to improve performance
- ✔ Automatic clog removal, intelligent flush & high voltage cleaning technologies
- ✔ High precision blood collecting system
- ✔ Fully automatic sample dilution in pre-diluent mode
- ✔ Fast and easy keyboard operation
- ✔ Cyanide free reagent to avoid environment contamination



LAB Market

لاب ماركت

- مستلزمات معملية
- كيمائيات
- أجهزة
- صيانة

01210312346
01016163354

بنها - عمارة ٢ غرب الاستاد - بجوار كنتاكي - خلف نقابة الاطباء ومطعم الرشيدى