

Emperor

EMP-168 BIOCHEMISTRY ANALYZER



☀️ FRIENDLY INTERFACE

- Windows CE. NET operation system
- Angle adjustable LCD Display and touch screen(optional)
- Inbuilt thermal printer and external printer(optional)
- USB mouse, virtual screen keyboard and inbuilt membrane keyboard
- Connect PC via RS-232 to transfer and print the data



☀️ MEASURE SYSTEM

- End point, kinetics, two points and absorbance
- Immunity nephelometry measure method
- Reagent blank and sample blank
- Multi-standards (linear or non-linear)
- Memory for reagent blank

☀️ FLOW CELL/CUVETTE

- 32µl quartz flow cell
- Inbuilt bottom position localizer
- Two-holder for cuvette/flow cell
- Peltier temperature control



لاب ماركت

- مستلزمات معملية
- كيمائيات
- أجهزة
- صيانة



LAB Market


01210312346
01016163354

بناها - عمارة ٢ غرب الاستاد - بجوار كنتاكي - خلف نقابة الاطباء - ومطعم الرشيدى



OPTICAL SYSTEM

- Halogen-tungsten lamp
- Automatic by 8 positions filter-wheel
- 7 filters(340nm~670nm) plus 1 free position



POWERFUL SOFTWARE

- Power on self-check and malfunction alarm
- Lamp automatic protection, user-define setting time
- Real time absorbance curve and temperature display
- Open system for reagent
- Support on-line printing
- Multiform reports by project, patient or history data
- Support user-define template test



لاب ماركت

- مستلزمات معملية
- كيمائيات
- أجهزة
- صيانة

LAB Market

01210312346
01016163354

بنها - عمارة ٢ غرب الاستاد - بجوار كنتاكي - خلف نقابة الاطباء ومطعم الرشيدى