

Global Diagnostics Supplier



## URIT-810

### Chemistry Analyzer

- High performance
- Easy operation & less maintenance



## لاب ماركت

- مستلزمات معملية
- كيمائيات
- أجهزة
- صيانة

# LAB Market

01210312346  
01016163354

بنها - عمارة ٢ غرب الاستاد - بجوار كنتاكي - خلف نقابة الاطباء ومطعم الرشيدى

# URIT

# URIT-810

## Chemistry Analyzer



### | Specifications |

**Analysis method:** Kinetic, end point, non-linear, sample blank, reagent blank, etc

**Light source:** halogen lamp, 6V/10W

**Wavelength (in nm):** 340,405,492,510,546,578,630

**Absorbance range:** -0.3~3.0Abs

**Wavelength accuracy:**  $\pm 1$ nm

**Spectrum band:**  $\pm 6$ nm

**Flowcell:** 10mm quartz cell, 32iL volume

**Temperature control:** Room temperature, 25°C, 30°C, 37°C, ( $\pm 0.1^\circ\text{C}$  variable)

**Sampling volume:** 100iL ~ 9999iL settable, recommended: 500iL

**Cross contamination:**  $\leq 1.0\%$

**Stability:**  $\leq 0.005$ Abs/h

**Memory:** Up to 320 items and 20000 test results

**Display:** Large LCD display

**Printer:** Built-in thermal printer

**Interface:** Serial port (RS232), Parallel port (For external printer)

**Ambient:** Room temperature: 10 °C ~ 32 °C; Humidity:  $\leq 85\%$

**Power:** AC (110-240)V, (50-60)Hz

**Dimension:** 392 x 375x 205mm ( length \* width \* height )

### | Features |

- High-quality filter & close optical system
- Convenient maintenance design
- Large memory capacity
- Advanced unexpected power-off protection

## URIT MEDICAL ELECTRONIC (GROUP) CO., LTD.

Head Office

Add.: No.4 East Alley, Jiuhua Road, Guilin,  
Guangxi 541001, P.R. China

Tel: +86-773-2288586 Fax: +86-773-2288560

E-mail: export@uritest.com service@uritest.com

Shenzhen Office

Add.: Room EJ, 24th Floor, New Baohui Building, Nanshan Road,  
Nanshan District, Shenzhen, P.R. China

Tel: +86-755-26050853 Fax: +86-755-86036150

Http://www.urit.com



CNAB047-Q