

AGD2020

Clinical Chemistry Analyzer



- Touchscreen Display
- Online Graphics
- Collated Patient Reports
- L - J Chart for Quality Control

Technical Specifications

● Assay Types

Absorbance, End point, Kinetic,
Fixed Time and Coagulation

● Aspiration Volume

From 0.3 mL to 1.0 mL by
specially designed Peristaltic Pump

● Printer

Built-in Thermal Printer which can be
switched ON/OFF, when desired.

● Wavelengths

Interference Filters 340, 405, 450, 505,
546, 578, 630 and 670 nm

● Cuvette Volume

500 μ L

● Interface

External Keyboard / Mouse

● Light Source

6V/10W, Tungsten Halogen lamp

● Storage Capacity

196 programs, 999 test records.

● Power Supply

230 V AC \pm 10%, 50Hz.

● Flow Cell

18 μ L

● Storage Condition

Temperature: 10° C to 60° C
Relative Humidity: up to 85%

● External Dimensions

360 mm (L) x 355mm (B) x 170mm (H)

● Temperature Control

25, 30 and 37° C by Peltier Effect

● Display

7" Graphic LCD with resistive
touchscreen,
Resolution - 800 x 480 pixels

● Weight

6.5 Kg (approx)



LAB Market

بنها - عمارة ٢ غرب الاستاد - بجوار كنتاكي - خلف نقابة الاطباء ومطعم الرشيدى

لاب ماركت

- مستلزمات معملية
- كيمووبات
- أجهزة
- صيانة

01016163354

CE

