

Wondfo®

Blood Gas Analyzer

Test Items

Basic blood-gas	pH, pCO ₂ , pO ₂
Haematology	Hct
Electrolyte	K ⁺ , Na ⁺ , Ca ⁺⁺ , Cl ⁻
Biochemistry	Glu*, Lac*

Test Packages

Test Card Parameters	BGDR-3	BGDR-8	BGDR7-Ca	BGDR7-Cl	BGDR6	BGDR3	BGDR4	BGDR10*
pH	*	*	*	*	*			*
pCO ₂	*	*	*	*	*			*
pO ₂	*	*	*	*	*			*
K ⁺		*	*	*	*	*	*	*
Na ⁺		*	*	*	*	*	*	*
Ca ⁺⁺		*		*			*	*
Cl ⁻		*	*		*	*	*	*
Hct		*	*	*				*
Glu*								*
Lac*								*



Wondfo®



BGA-101

Blood Gas Analyzer

- Quickly
- Portable
- Accurate
- Efficient

GUANGZHOU WONDFO BIOTECH CO., LTD.



LAB Market

لاب ماركت

- مستلزمات معملية
- كيمائيات
- أجهزة
- صيانة

01096411333

بنها - عمارة ٢ غرب الاستاد - بجوار كنتاكي - خلف نقابة الاطباء ومطعم الرشيدى

Wondfo®


Blood Gas Analyzer






Wondfo®



Operation



Contents



لاب ماركت

- مستلزمات معملية
- كيمائيات
- أجهزة
- صيانة

LAB Market

01096411333

بنها - عمارة ٢ غرب الاستاد - بجوار كنتاكي - خلف نقابة الاطباء ومطعم الرشيدى