



**AADEE** S.A.

**RUMI**  
Blood Gas and pH Analyzer



**Measured Parameters**

pH, pCO<sub>2</sub> y pO<sub>2</sub>.

**Calculated Parameters**

HCO<sub>3</sub><sup>-</sup>, ST.pH, Bstd, BEp, EBs,  
BEfec, O<sub>2</sub>sat, TCO<sub>2</sub> y O<sub>2</sub>CT

RUMI BG has a highly versatile patented device that allows 22 /30-per-hour samples input by syringe or capillary aspiration, or syringe injection.

