

لاب ماركت

- مستلزمات معملية
- كيمويات
- أجهزة
- صيانة

icon

LAB Market

01096411333

بنها - عمارة ٢ غرب الاستاد - بجوار كنتاكي - خلف نقابة الاطباء ومطعم الرشيدى

A Smart Hematology Analyzer

- Smallest sample volume makes it ideal for pediatric use
- Smallest reagent consumption
- A4 paper size footprint
- ER mode for use at the point-of-care
- NormaAccess – Remote support tool



TECHNICAL DATA

Parameters

WBC, LYM, MID, GRA, LYM%, MID%, GRA%,
RBC, HGB, MCV, HCT, MCH, MCHC, RDWcv, RDWsd,
PLT, MPV, PCT, PDWcw, PDWsd, P-LCR%, P-LCC

Measurement technology

Volumetrical impedance method for WBC and RBC
Photometrical method for HGB

Histograms

WBC, RBC, PLT

Sampling mode

Closed and open vial mode

Aspirated sample volume

Closed vial: 38 µl

Open vial: 13 µl

Processed sample volume

2.5 µl

Throughput

60 tests per hour

Precision (CV%)

WBC < 3%, RBC < 2%, PLT < 5% MCV < 1%, HGB < 2%

Storage capacity

40,000 test results incl. histograms

User interface

LCD, 10.1", 1280 x 800 with capacitive touchscreen, portrait

Connectivity

2 x USB ports, Ethernet, LIS (HL7 and ASTM)

Dimensions (HxWxD)

280 x 216 x 320 mm

Net weight

8.7 kg

Norma Icon is ideal for doctors' offices, emergency departments, small to medium size labs and polyclinics.



- Microfluidic technology
- Patented, innovative solutions
- Trendsetting design



Shortlisted for the Senetics Innovation Award 2014
Category "Most Innovative Product in the Health Care Sector"

reddot award 2014
winner



icon